

Report Qualisys Motion Capture on ROS and Network Configuration

Created by Antonio Camposeo - Polytechnic of Bari

January 2024

1 Qualisys Motion Capture on ROS

1.1 Virtual Machine Installation

First, make sure that you have installed the Qualisys software on your PC and have activated the licence, then it is possible to follow this step to get the data given by the software on Windows directly on Linux/Ubuntu.

Once you have checked that the software is working correctly, make sure that the cameras are calibrated and that the programme knows the position of the body to be tracked. We will now proceed with the installation of Ubuntu LST 20.04 on Windows, followed by the installation of ROS on the virtual machine.

Step for Windows 10 (*Check if these steps are the same for Windows 11*):

- Open Microsoft Store;
- Search Ubuntu 20.04.06 LTS on the search bar and download it.




Figure 1: Install Ubuntu 20.04.6 selected in the red box

- Wait the installation and open it.

After the correct installation of Ubuntu20.04.6 on Windows10, open the installed programme and follow the procedure to install Ubuntu on it. At the end you will be asked to enter your username and password. If there are any problems, go to the section Possible Problems.

1.2 Possible Problems

One of the possible problems could be the error shown in the picture below.




Figure 2: Ubuntu Installation Problem

To fix this problem, go on the search bar of windows and find Windows Features. Open it and check if the box "Virtual Machine Platform" is checked, if not, select it in. Then restart the computer and open Ubuntu20.04.6 LTS.

The installation will start correctly, and then you will be able to use the virtual machine on Windows.

1.3 Ros Installation on Virtual Machine

To install Ros Noetic on Ubuntu20.04.6 LTS, follow the command list below.

Commands list

#Package Setup

```
sudo sh -c 'echo "deb http://packages.ros.org/ros/ubuntu $(lsb_release -sc) main"
> /etc/apt/sources.list.d/ros-latest.list'
```

#input ROS key

```
sudo apt install curl # if you haven't already installed curl
curl -s https://raw.githubusercontent.com/ros/rosdistro/master/ros.asc
| sudo apt-key add -
```

#Update

```
sudo apt update
```

#Install

```
sudo apt install ros-noetic-desktop-full
```

#Source environment

```
source /opt/ros/noetic/setup.bash
```

1.4 Qualisys Motion Capture SDK on ROS

Once the ros noetic installation is done, it is possible to create a ros package.

#Create project directory

```
mkdir QROS2
```

```
#Move to project directory  
cd QROS2
```

```
#Create source directoty within project directory copying the Qualisys source driver  
mkdir src
```

```
cd src
```

```
#clone example project to source repository  
git clone https://github.com/KTH-SML/motion_capture_system.git
```

```
#return to project folder  
cd ..
```

```
#build driver  
catkin_make
```

```
#source new build environment  
source devel/setup.bash
```

Now, test if the project insert work.

```
#In the same terminal  
#launch Qualisys ROS node  
roslaunch mocap_qualisys qualisys.launch server_address:=192.168.56.1 udp_port:=-1
```

```
#open new window terminal  
#navigate to directory folder  
cd QROS  
#souce build environment  
source devel/setup.bash  
#check if there are topics  
rostopic list
```

1.4.1 Useful commands ROS

```
#tools  
rostopic #shows rostopic commands
```

```
rostopic list #shows all published nodes
```

```
rostopic info #shows info about node
```

```
rostopic echo #shows data from node
```

```
#new window to pull data from any location  
#navigate to directory folder  
cd QROS
```

```
#source build environment  
source devel/setup.bash  
#show what is using each node  
rostopic info /qualisys/name_of_tracking_obj/pose #example
```

```
#plot data  
export DISPLAY=$(ip route list default | awk '{print $3}'):0 #opens new window
```

```
export LIBGL_ALWAYS_INDIRECT=1 #creates plot
```

```
rqt_plot /qualisys/name_of_tracking_obj/pose/pose/position/z #plot that looks at ground
```

2 Network Configuration

A ROS Network had to be used to interconnect all devices such as the Ground Control Station, the Rov, the haptic system and the motion capture system. Through the ROS environment, it was possible to communicate with any node in the network at any time. The requirements to set up such a network are:

- Complete two-way connectivity between all pairs of machines, on all ports.
- Each machine must publish itself with a name that all other machines can solve.

To establish full connectivity, it was necessary to create a local network, in which all devices are connected to an Ethernet switch. Each device is assigned a static IP to remain stable over time. Therefore, the devices in the network are configured as follows:

- Ground Control Station : IP 192.168.2.1
- Qualisys Motion Capture : IP 192.168.2.3
- BlueROV2 Fathom X : IP 192.168.2.2




Figure 3: System Local Network

After configuring the Ros Network, we must define the Ros Master. For simplicity and efficiency, we use the Control Station as the master node, as it allows us to process data faster and have control over the entire network.

2.1 Static IP Configuration

As mentioned at the beginning of this section, each device connected to the network must have a static IP. Let us now set the ip for each device, knowing that the default ip of the BlueROV2 system is 192.168.2.2.

2.1.1 Motion Capture Computer - Operating System : Windows10

2.1.2 Control Station Computer - Operating System : Ubuntu20.04

2.1.3 Check interconnection between devices

<http://wiki.ros.org/ROS/NetworkSetup>

To check that the connection between devices is working, run on terminal of each devices this command:

```
ping 192.168.X.X
```

192.168.X.X is the static ip related to the other device in the net. If the data are transmitted and received the device is connected.

Check the connection in both devices.

```
# from ground control station
ping 192.168.2.3
# from qualisys computer
ping 192.169.2.1
```

2.2 ROS Network Configuration

Ros Master:

```
export ROS_IP=ip_of_this_machine
export ROS_MASTER_URI=http://ip_of_this_machine:11311
```

Ros Node:

```
export ROS_IP=ip_of_this_machine
export ROS_MASTER_URI=http://ip_of_master:11311
```

2.3 OpenCM9.04 Access via Ground Control Computer

In order to access the opencm card using the Ethernet connection of the bluerov system, third party software installed on the RPi3 was used. The software is called VirtualHere, and it allows us to virtualise the serial connection between the RPi3 and a USB device, giving us the IP address to access that device. So we can control it without it being physically connected.




Figure 4: VirtualHere Diagram

RPi3

Download <https://www.virtualhere.com/sites/default/files/usbserver/vhusbdarm> and install it on RPi3, and write a bash script in order to run it every time the operating system starts.

Ground Control Computer : Ubuntu20.04

Download <https://www.virtualhere.com/sites/default/files/usbclient/vhuit64> and install it on ubuntu using these commands:

```
# Go to the folder containing the file  
sudo chmod +x ./vhuit64  
# Run it  
sudo ./vhuit64
```